

"Quality has been proven 50 countries'







Clean room facility



Tensile strengthmeasurement



Moisture measuring system



Cog machine



Special cog machine



-1-



NVL-Bio lifting threads are made of PDO, PCL which are bio-degradable materials. This

and promote skin regeneration.

NVL-Bio Lifting threads have various types and sizes for different applications.

[Mono/ Twin/ Screw/ Double Screw/ Hydra/ Coil spring/ Braid/ Cog/ Double Cog/ Nose/

Eye/ Body/ Pro]

Especially our Double cog threads are designed by patented technology for face-lifting applications. Unlike other threads whose surfaces are cut or knotted or just mold to create cogs, the cogs of NVL_Bio Double Cog threads are created by mold-cutting techniques without heating the threads. This technology minimizes the damage to the surfaces of the threads and makes the threads solid therefore, its tensile and anchoring strengths are superior to other thread products by patented unique technology.

TYPE	PDO	PCL	PLCL
Biodegradation Period	6-7 months	36 months	18-23 months
V	107		
Early formation	After	30 days	After 90 days
		2-	

NVL-BIO lifting THREAD





MONO

NEEDLE (GUAGE)	NEEDLE (mm)	USP
26G	38	5-0
26G	60	5-0
26G	90	5-0
27G	38	5-0
27G	50	5-0
29G	38	6-0
29G	50	6-0
29G	60	6-0
30G	25	7-0
30G	38	7-0





TWIN

NEEDLE (GUAGE)	NEEDLE (mm)	USP	
26G	38	5-0	
26G	60	5-0	
269	90	5-0	

-4-

NVL-BIO lifting THREAD







SCREW

NEEDLE (GUAGE)	NEEDLE (mm)	USP
26G	38	6-0
26G	60	6-0
27G	38	6-0
27G	50	6-0
29G	38	6-0
29G	50	6-0





DOUBLE SCREW

NEEDLE (GUAGE)	NEEDLE (mm)	USP	5
26G	38	7-0	
26G	60	7-0	
26G	90	7-0	
27G	38	7-0	
27G	50	7-0	
27G	60	7-0	











PDO


NEEDLE	NEEDLE	SUTURE	
(GUAGE)	(mm)	(mm)	A comment
266	38	30	
29G	50	40	





DOUBLE COGI (MOLDING PRESSES COGS)






NVL-BIO lifting THREAD






DOUBLE COGII (MOLDING PRESSES COGS)








NVL-BIO lifting THREAD


100

12D





2-0



PRO (DOUBLE NEEDLE)

N. O	PDO	PGL					
NEEDLE (GUAGE)	NEEDLE (mm)	SUTURE (mm)	USP	NEEDLE (GUAGE)	NEEDLE (mm)	SUTURE (mm)	USP
196	120	30	2	196	100	850	2
19G	120	30	2	19G	120	280	2
216	110	290	1-0	196	120	450	2
21G	110	440	1-0	19G	120	850	2

-12-



- LIFTING POWER -



K DERMA CO., LTD.

#201 , 29-5 Nonhyeon-ro95-gil, Gangnam-gu, seoul, south korea Tel: +82-10-493-0891 e-mail: info@kderma.com WWW.KDERMA.COM